


The use of SomnoFit in case of overbite (mandibular retrognathia)

For those who suffer from severe overbite, the standard bands of a snoring splint may be too short in individual cases. This is because the bands are used to adjust the advancement of the lower jaw; if the lower jaw is naturally set very far back, even the longest band length would pull it so far forward that this could be uncomfortable to wear.



With SomnoFit, this problem can be elegantly solved. Take band no. 6 (= the longest band) and cut it into two parts as shown in Figure 2. This way you get two separate bands, which you then hook into the SomnoFit as usual (Figure 3). These bands are extra long and cause the lower jaw to be pulled forward only very moderately. This is the ideal solution for people with a severe overbite.

